THE PIONEER OF MEDICAL EDUCATIONAL SIMULATOR +

Airway Management Simulator BT-CSIE

The airway management simulator with an objective feedback system

NURSING MEDICINE MILITARY EMS HOSPITAL





Airway Management Simulator

MANIKIN

- · Lifelike and durable silicon skin
- · Head tilt-Chin lift Maneuver, Jaw thrust maneuver, Sniffing Position
- NPA, OPA, ETT, LMA/i-gel, Combi-tube
- Visible Chest rising / Gastric inflation during ventilation
- Breathing sound with stethoscope during ventilation

SOFTWARE

- Changeable Guideline
- 3 MODES (ETT, LMA/i-gel, Combi-tube)
- · Alarms for incisors due to laryngoscope
- Displaying the location of the tube (ETT)
 - · Tracheal intubation (Correct/Excessive), Gastric intubation
- Real-time feedback such as volume, speed etc. during ventilation
- Objective evaluation with auto/manual checklists: Scoring
- Saving and Printing out results



* Training screen



COMPONENTS

Adult Torso: 1EATablet: 1EA

Power Adaptor: 1EABattery holder: 1EAStorage case: 1EA

SPECIFICATION

Dimensions: 70 x 30 x 22cm

Weight: 10kg



