

COMMODITY SPECIFICATION

Product	CPR Evaluation Simulator	Model	BT-CPTA PLUS
Manufacture	BT INC	Q'ty	1 SET



A. FEATURES

- Human rib-based spring, Lifelike resistance.
- Articular arms and legs ; possible to change the position.
- Wireless connection between simulator and operating device; Possible to evaluate up to 6 simulators at once.
- Displaying the multi-screen while operating multiple simulators.
- Detect up to 7cm of compression depth ; Can guide the trainees for the precise chest compression.
- Offering the magnetic sticker for AED pads ; Possible for the semi-permanent training of AED pads.
- Set your own guideline, evaluation list, check list and etc. for the evaluation.
- After the training, detailed result and the analysis data can be saved/printed out; these can be exported to Excel.
- Possible to check the airway opening, compression depth, hand position, compression rate, the number of compressions, hands-off time, ventilation volume, the number of ventilations, expiratory time, expiratory rate and the accuracy in real time.
- Displaying the palpation time when the carotid artery is palpated. Display each hands-off time separately on every hands-off section.
- Providing compression fraction data in result to check the actual the time for the compression during CPR.
- Metronome function provides 100, 110, 120 bmp and possible to turn on/off the carotid artery according to the scenario.
- For the power, AC 220V is used to operate the manikin. In case the power is cut, the result data would be

COMMODITY SPECIFICATION

preserved by using the battery.

B. SPECIFICATION

► Display windows

- Trainee information : Save the group&trainee information (Possible to set the number of trainee in group)
- Setup for the training: Set up the guideline, scenario, evaluation list and points, training time and checklist.
- Training screen :
 - Displaying trainee info. / Training time / Status of manikin / Chest compression and ventilation data
 - Check consciousness : After checking the patient's consciousness, it will be marked at graph & history.
 - Check pulse : After checking the patient's carotid pulse, it will be marked at graph & history.
 - Chest compression : Displaying the detailed data such as compression depth, hand position, incomplete recoil, compression rate, hands-off time and the number of compression in graph and accuracy analysis during the chest compression.
 - Airway opening : When the airway is opened, it will be marked at graph & history. (Also it detects Hyper-extension)
 - Artificial respiration : Displaying the detailed data such as ventilation volume, expiratory time, expiratory rate and the number of ventilation at graph and accuracy analysis during the ventilation.
 - Defibrillation : When virtual shock is delivered, shock delivered time, the number of shock will be marked in graph & history / Information of ECG, pulse, RR will be displayed on the manikin's status.
 - Virtual patient monitor : Showing the patient's vital signs when there is an extended monitor.
- Accuracy analysis : Displaying analysis data for chest compression, ventilation as a chart in real-time.
- Training record : All practice records from beginning to the end will be recorded in graph and history.
- Checklist : The preset checklist for the evaluation will be displayed.
- History : Showing the history when auto and manual checklist items are marked.
- Metronome : Set up the metronome. (OFF, 100bpm, 110bpm, 120bpm)
- Scenario : Displaying ECG, pulse, RR according to the preset information.
- Evaluation result :
 - Displaying Trainee info and result data
 - Compression : The number of total compression, average depth, rate, compression fraction, compression per minute, hands-off-time(separately displayed for each hands-off), the number of each good/insufficient/excessive compression, the number of each good/upper/lower/right/left compression, the number of complete/incomplete recoil.
 - * All data is converted to a percentage(%) and displayed with the information
 - Ventilation : The number of total ventilation, ventilation per minute, average volume,

COMMODITY SPECIFICATION

average expiratory rate,

the number of each good/insufficient/excessive expiratory rate and gastric ventilation,

the number of each good/insufficient/excessive ventilation volume,

the number of each good/insufficient/excessive expiratory time.

- Displaying the accuracy of each evaluation list, score and the total score.
- Displaying the time of each checklist done, the score and the total score.
- Save/Print the evaluation result

– Previous result :

- Display the each trainee's previous result.
- Result summary / Training screen/ Accuracy analysis / History's selective display.
- Print/Export the result data in Excel format.

C. CONSIST OF

- CPR Training manikin (Adult Torso) – 1EA
- Power adaptor (12VDC 1.5A) - 1EA
- Battery Holder - 1EA
- Articulated arms and legs – 1 SET
- Magnetic stickers for AED pads – 1 SET
- Storage case – 1 SET
- Operating Software included

D. WEIGHT & SIZE

- CPR training manikin(Whole body) : 18kg (168×50×21cm)
- Storage case : 9.2kg (90×40×40cm)
- Total : Approximately 29kg