

COMMODITY SPECIFICATION

Product	Male Catheterization Simulator	Model	BT-CSCM
Manufacture	BT INC	Q'ty	1 SET



A. FEATURES

- Life-like practical training
- Make position similar to the real patient
- Any water-soluble Jelly being used at clinical situation can be applied
- Training with various size of catheter is available
- Non-return valve : Prevent leaking after taking catheter out.
- Checking fluid after correct insertion of catheter at the right position.
- Foreskin (attachable / detachable)
- Genital organs and buttocks are attachable : easy to maintain

B. CONSIST OF

- Adult buttocks model : 1set
- Simulated water bag : 1EA
- Storage case : 1EA

C. WEIGHT & SIZE

- Adult buttocks model : 3.5kg (430×280×200mm)
- Storage case : 2.4kg (440×330×390mm)
- Total : Approximately 6kg