

COMMODITY SPECIFICATION

Product	Enema Administration Simulator	Model	BT-CSEN
Manufacture	BT INC	Q'ty	1 SET



A. FEATURES

- Life-like practical training with BT Product
- Make the Sim's Position, Knee-chest Position
- Any water-soluble Jelly being used at clinical situation can be applied
- Any Glycerin being used at clinical situation can be applied
- Training with various size of catheter is available
- Non-return valve : Prevent leaking after taking catheter out
- Artificial feces can be inserted / removed

B. CONSIST OF

- Adult Buttocks model : 1set
- Drain tube : 1EA
- Storage bag : 1EA

C. WEIGHT & SIZE

- Adult Buttocks model : 3.5kg (330×430×340mm)
- Storage bag : 2.4kg (440×330×390mm)
- Total : Approximately 6kg