

# COMMODITY SPECIFICATION

<b>Product</b>	<b>Intramuscular Injection Trainer</b>	<b>Model</b>	<b>BT-CSIM-T</b>
<b>Manufacture</b>	<b>BT INC</b>	<b>Q'ty</b>	<b>1 SET</b>



## A. FEATURES

- The touch feeling of the model is similar to human skin.
- The buttocks model is composed of epidermis layer, muscular layer and skeletal frame.
- Ventrogluteal and dorsogluteal site injection practice; both side
- The skin is made of soft silicon that leaves no needle marks with self-healing system
- Both side of buttocks can be used for practice.
- Fluid injection and drain
- LED lights / Warning Sound can show the feedback in real time
- Power : AC power or battery

# COMMODITY SPECIFICATION

## B. CONSIST OF

- Adult buttocks model - 1EA
- Buttocks holder - 1EA
- Power Adaptor (9V) - 1EA
- Battery Holder - 1EA
- Drain bottle (including one drain tube) - 1EA
- Drain Syringe - 1EA
- Storage Case - 1EA

## C. WEIGHT & SIZE

- Adult buttocks model: 5.5kg (330×260×380mm)
- Storage Case : 2.7kg (450×450×260mm)
- Total : Approximately 8kg