

# COMMODITY SPECIFICATION

<b>Product</b>	<b>Intramuscular Injection Simulator</b>	<b>Model</b>	<b>BT-CSIM1</b>
<b>Manufacture</b>	<b>BT INC</b>	<b>Q'ty</b>	<b>1 SET</b>



## A. FEATURES

- The touch feeling of the model is similar to human skin.
- The buttocks model is composed of epidermis layer, muscular layer and skeletal frame.
- Ventrogluteal and dorsogluteal site injection practice; both side
- The skin is made of soft silicon that leaves no needle marks with self-healing system
- Fluid injection and Automated drain
- Advanced Education method with touch screen & color display.
- Self-Training and Evaluation of user's performance with Training & Evaluation mode.
- Automated on-line connection with Bluetooth interface.
- Objective assessment for procedures before Injection and injection
- All results being displayed in real time
- The detailed evaluation results are stored in JPEG format
- Power : AC power or battery

# COMMODITY SPECIFICATION

## B. SPECIFICATION & CONSIST OF

- Adult buttocks model - 1EA
- Tablet PC - 1EA
- Buttocks holder - 1EA
- Drain bottle (including one drain tube) - 1EA
- Power Adaptor (9V) - 1EA
- Battery Holder - 1EA
- Storage Case - 1EA

## C. WEIGHT & SIZE

- Adult buttocks model: 5.5kg (330×260×380mm)
- Storage Case : 2.7kg (450×450×260mm)
- Display device : Variable depending on situation
- Total : Approximately 10kg