
# **COMMODITY SPECIFICATION**

Product	Intravenous Injection Training Infants Arm Model	Model	BT-CSIV-I
Manufacture	BT INC	Q'ty	1 SET



## A. FEATURES

- IV injection, blood collection, IV Route Status Practice is possible at 2 points of back hand and 1 point of popliteal
- Easy to attach and detach arms to the arm mount with magnetic holder
- Consist of 2 types of arms; 1-year-old, 3-years-old
- Flexible finger and easy to change the angle of arm.
- Can check flash-Back technique
- Hybrid system; Circulation pump and gravity
- Circulation pump: easy to set up, clean and maintain
- Power : AC power or battery



# **COMMODITY SPECIFICATION**

## **B. CONSIST OF**

- Infant Arm Model(Vein Tube installed) 2EA
- Motor Pump 1EA
- Timer Switch 1EA
- Arm Holder Adapter 1EA
- Mount 1EA
- Simulated Blood (30ml) 1EA, Plastic Container 1EA, Connection Tube 3EA
- AC adapter 1EA
- Battery Holder 1EA
- Storage case 1EA

### C. WEIGHT & SIZE

- Arm model 1 (3 years old): 0.7kg (370×100×80mm)
- Arm model 2 (1 years old) : 0.5kg (250×80×80mm)
- Circulation pump :  $0.6 \mathrm{kg}$  (180×70×105 mm)
- Mount : 0.3 kg (270×200×24mm)
- Storage Case : 2.2kg (730×280×140mm)
- Total : Approximately  $5.5 \mathrm{kg}$

