

COMMODITY SPECIFICATION

Product	Intravenous Injection Training Model	Model	BT-CSIV-S
Manufacture	BT INC	Q'ty	1 SET



A. FEATURES

- Lifelike and realistic arm model
- Soft skin silicon ; less needle marks after training.
- Injection and blood draw practice on Basilic vein, Cephalic vein, Median cubital vein
- One touch stand ; easy to set up
- Checking flashback

B. CONSIST OF

- Arm Model - 1EA
- Stand – 1EA
- Simulated Blood – 1EA
- Blood supply bottle & Rubber stoppers – 2EA
- Storage Case – 1EA

C. WEIGHT & SIZE

- Arm Model : 0.2(230×90×65mm)
- Stand : 0.5kg (350×100×50mm)
- Total : Approximately 1kg