

COMMODITY SPECIFICATION

Product	Intravenous Injection Training Arm Model 1	Model	BT-CSIV1
Manufacture	BT INC	Q'ty	1 SET



A. FEATURES

- Skin pad is consisted of epidermis, subcutaneous fat, vein.
- Injection practice on basilic vein, cephalic vein, median cubital vein and metacarpal veins
- All-in-One skin pads(skin and vein) ; easy to replace
- One touch connecting method ; easy to set up
- Hybrid system ; Circulation pump and gravity
- Circulation pump : easy to set up, clean and maintain
- Power : AC power or battery

COMMODITY SPECIFICATION

B. CONSIST OF

- Left Arm Manikin - 1EA
- Skin Pad (Vein tube installed) - 2EA
- Motor Pump - 1EA
- Arm Holder – 1EA
- Simulated Blood(30ml) - 1EA, Blood Container – 1EA, Connecting Tube - 3EA
- Power Adaptor (9VDC 1A) - 1EA, Battery Holder - 1EA
- Storage Case – 1EA

C. WEIGHT & SIZE

- Left Arm Model : 1.4kg (580×120×80mm)
- Circulation pump : 0.6kg (180×70×105mm)
- Arm Holder : 0.3kg (270×200×24mm)
- Storage Case : 2.2kg (730×280×140mm)
- Total : Approximately 5.5kg