

COMMODITY SPECIFICATION

Product	Arterial Puncture Training Model	Model	BT-CSPA
Manufacture	BT INC	Q'ty	1 SET



A. FEATURES

- Arm model and skin pad is similar to real human arm
- Palpate life like pulse at radial artery
- Changeable pulse intensity and rate
- Check the backward flowing of blood while arterial puncture training
- Realistic resistance at the arterial wall during arterial puncture training
- Hard to find needle marks after punctures
- Skin&Vein All-in-one pad; Easy installation and replacement of skin pad
- The color of palm is changing while Allen test
- Easy clean, Easy maintenance; Self-cleaning system
- User can select AC power or Battery mode

COMMODITY SPECIFICATION

B. CONSIST OF

- Adult arm manikin
 - Arm manikin : 1ea
 - Skin Pad : 2ea
- Blood circulation system : 1ea
- Blood supply bottle : 1ea
- Simulated blood : 1ea
- Power adapter : 1ea
- Storage case : 1ea
- User manual : 1ea

C. WEIGHT & SIZE

- Arm Model : Approximately 0.6kg (110×265×95mm)
- Blood supply system : approximately 2kg (307×307×170mm)
- Total model : Approximately 3kg