

COMMODITY SPECIFICATION

Product	CPR&Airway Management Training Simulator	Model	BT-SEEM-AIR
Manufacture	BT INC / South Korea	Q'ty	1 SET



A. FEATURES

- Training Chest Compression, Airway management, Ventilation for basic life support
- Intubation by using laryngoscope
- Insert the oral airway (OPA) and nasal airway(NPA) to the manikin
- By attaching offering magnetic stickers to AED trainer pads, AED trainer pads can be used semi-permanently
- Feedback about chest compression is being displayed through LED panel at the abdomen
- Training and evaluation of CPR with S/W
- Available to check evaluation result and print it out with S/W
- Self-training of CPR with a tablet or laptop / Group training with up to 30 manikins at the same time

B. SPECIFICATION

■ Airway management

- Head tilt -Chin lift, Jaw thrust maneuver
- Ventilation with a bag-valve-mask
- Insert OPA/NPA
- Intubation with a laryngoscope
- Intubation with an endotracheal tube (ETT), laryngeal mask (LMA/I-gel), combi tube (Combi-Tube)
- Observable chest rising during the ventilation
- Click sound for the front teeth due to the excessive use of laryngoscope

COMMODITY SPECIFICATION

■ Chest compression

* Wireless connection between manikin and tablet/smartphone via Bluetooth

- Easy-to-use and user-friendly S/W
- Possible to evaluate up to 30 manikins at the same time ; tablet, printer and laptop are optional
- With its S/W, change the guide line before training
- With its S/W, display the compression depth, compression rate, hand position and also the ventilation volume
- With its S/W, check the result after training and save/print it
- Export the result data (Result in statistics) to excel format
- Score the training result and rank the saved scores
- Choose the mode between Training mode and Evaluation mode

* LED Display embedded in the manikin

- Real-time feedback of chest compression depth and rate with LED lamp
- Display total number of compression and the number of good compression
- Display the instant depth of compression
- Detect incomplete recoil
- Detect the good chest compression when depth and rate are good and give sound feedback (Sound On/Off)
- Display the feedback about ventilation such as volume and expiratory time
- Display the total number of ventilation and the number of good ventilation
- Display the ventilation volume

* Manikin body

- Observable chest rising during ventilation
- By attaching offering magnetic stickers to AED trainer pads, AED trainer pads can be used semi-permanently
- The storage bag can be used as a training mattress
- Use batteries / AC adapter

COMMODITY SPECIFICATION

B. COMPONENTS

- Adult Torso manikin : 1EA
- Storage case : 1EA
- Lubricant for intubation : 1EA
- AC Power Adaptor, User Manual : 1EA

C. WEIGHT & SIZE

- Adult Torso : 6.4kg (670×330×210mm)