

# COMMODITY SPECIFICATION

<b>Product</b>	<b>Self-training CPR Model</b>	<b>Model</b>	<b>BT-SEEM-P</b>	
<b>Manufacture</b>	<b>BT INC</b>	<b>Q'ty</b>	<b>1</b>	<b>SET</b>



## A. FEATURES

- ✓ **Check the feedback with the display indicator built into the main body**
  - Real-time confirmation of chest compression depth and speed with LED lamps
  - Can check the total number of chest compressions and the exact number of chest compressions
  - Can display chest compression depth values
  - Possible to display chest compression position (up/down/left/right)
  - Incomplete recoil can be checked
  - Good chest compressions (both compression depth and normal speed) can be confirmed by sound
  - Real-time confirmation of ventilation volume and expiratory time with LED lamps
  - Can check the total number of ventilations and the exact number of ventilations
  - Can display ventilation volume
  - Good ventilation (both volume and expiratory time) can be confirmed by sound
  - Possible to ON/OFF display feedback & sound feedback

# COMMODITY SPECIFICATION

## ✓ Main body

- Possible to check the visual chest rising during ventilation
- Training for AED trainer semi-permanently with included self-adhesive rubber steel sheets
- Hygienic training is possible because the facial skin and lung bags are easy to be replaced
- The storage bag can be used as a training mattress
- Portable auxiliary battery (sold separately), batteries and AC power adapter can be used

## B. CONSIST OF

- Adult Torso : 1EA
- Storage case : 1EA
- Face skin: 1EA
- Lung bag: 2EA
- Self-adhesive Rubber steel sheet : 2EA
- AC Power Adaptor : 1EA
- User Manual : 1EA

## C. WEIGHT & SIZE

- Adult Torso : 4.4kg (650×335×215mm)
- Storage case : 1.5kg (670×350×230mm)
- Total : Approximately 6.2kg