

THE PIONEER OF MEDICAL  
EDUCATIONAL SIMULATOR +

# Intramuscular Injection Simulator

BT-CSIM1

Precise result according to the injection site and depth

NURSING MEDICINE MILITARY EMS HOSPITAL



# Intramuscular Injection Simulator

## MANIKIN

- Durable silicon skin & Easy to replace
- Fluid injection & Automatic drain with an embedded pump
- LED Landmarks for palpating bones
- Ventrogluteal and Dorsogluteal site (Left&Right)



Ventro/  
Dorso



Fluid  
Injection



LED  
Guide



Anatomical  
Skeleton




Positioning  
Holder




Self-  
changeable  
Skin

## SOFTWARE

- Tablet PC ► Evaluate practical training
- User friendly interface – Easy to use
- Training Mode & Evaluation Mode
- Detailed feedback in real-time
- Objective evaluation with a manual checklist : Scoring
- Saving and Printing out results
- Evaluate students according to the process



\*Display



\*Video guide & Score



\*Result

## COMPONENTS

- Buttocks body model : 1EA
- Tablet : 1EA
- Buttocks holder : 1EA
- Drain bottle: 1EA
- Drain tube : 1EA
- Power Adaptor : 1EA
- Storage case : 1EA

## SPECIFICATION

- Buttocks model : 33 x 26 x 38cm (5.5kg)
- Display device : 26 x 18 x 1cm (0.6kg)
- Total : approximately 10kg

