

THE PIONEER OF MEDICAL
EDUCATIONAL SIMULATOR +

Intravenous Injection Infants Arm Model

The advanced IV injection training model with a motor system

BT-CSIV-I

NURSING MEDICINE MILITARY EMS HOSPITAL



IV Injection Training Arm Model - Infants

MANIKIN

- Arm model is similar to Real Human Arm
 - Soft material with realistic tactile feeling
 - Practice IV injection, Draw blood, IV Route status at 2 points of back hand, 1 point popliteal
- Realistic resistance feeling on the skin and venous wall
- Easy Attach / Detachment and adjust angle easily
 - Flexible finger and easy to change the angle of arm
 - Easy to attach / detach with magnetic type of arm mount
- Pump-driven system
 - Battery mode – can be used wirelessly
 - Easy to set up
- Flash-Back
 - Blood flow is similar to human



Skin & Vein
All-in-one pad



Life-like
resistance
feeling



3 layers
(Epidermis, Vein,
Subcutaneous fat)



Durable
Skin



Motor



Flashback



Easy
Cleaning



Battery



Easy
Installation

COMPONENTS

- Baby arms (around 1-year-old) : 2EA
- Infant arms (around 3-year-old) : 2EA
- Circulation pump : 1EA, Stand : 1EA
- Blood supply bottle : 1EA
- Power Adaptor : 1EA
- Battery Holder : 1EA
- Storage case : 1EA

SPECIFICATION

- Arm model 1 : 37 x 10 x 8 cm (0.7kg)
- Arm model 2 : 25 x 8 x 8 cm (0.5kg)
- Mount : 2.7 x 2 x 2.4 cm (0.3kg)
- Circulation pump : 18 x 70 x 11 cm (0.6kg)
- Storage case : 73 x 28 x 14 cm (2.2kg)

