


THE PIONEER OF MEDICAL
EDUCATIONAL SIMULATOR +

Intravenous Injection Training Arm Model 1

The advanced IV injection training model with a motor system

BT-CSIV1


NURSING MEDICINE MILITARY EMS HOSPITAL



IV Injection Training Arm Model 1

FEATURES

- Lifelike arm model and skin pad
 - Training IV access, vascular ac
 - Practice on Basilic vein, Cephalic vein, Median cubital vein and Metacarpal vein
- Realistic resistance feeling on the skin and venous wall
 - Skin pad is composed of Epidermis / Subcutaneous fat / Vein



Skin & Vein
All-in-one pad



Life-like
resistance
feeling



3 layers
(Epidermis, Vein,
Subcutaneous fat)



Durable
Skin

- Easy installation & replacement of skin pad
 - All-in-one skin
 - Convenient to replace the skin pad
- Pump-driven system
 - Battery mode – can be used wirelessly
 - Easy to set up
- Flash-Back
 - Blood flow is similar to human



Motor



Flashback



Easy
Cleaning



Battery



Easy
Installation

COMPONENTS

- Arm Model : 1EA
- Skin pads: 1SET (2EA)
- Circulation pump : 1EA, Stand : 1EA
- Blood supply bottle : 1EA
- Power Adaptor : 1EA
- Battery Holder : 1EA
- Storage case : 1EA

SPECIFICATION

- Arm model : 58 x 12 x 8 cm (1.4kg)
- Circulation pump : 18 x 70 x 11 cm (0.6kg)
- Storage case : 73 x 28 x 14 cm (2.2kg)

