

THE PIONEER OF MEDICAL  
EDUCATIONAL SIMULATOR +

# Sherpa X, Smart CPR Training Model

CPR Training with accurate feedback in real-time and analysis

BT-SEEM2

NURSING MEDICINE MILITARY EMS HOSPITAL



# Sherpa X - Smart CPR Training Model

## MANIKIN

- LED Display for real-time feedback (Chest compression, Ventilation)
- Chest compression – Total number, Good number, Depth, Speed, Incomplete recoil, Hand position
- Ventilation – Total number, Good number, Volume, Rate
- Sound, LED ON / OFF
- Head tilt-Chin lift Maneuver, Airway opening
- Visible Chest rising during ventilation
- Semi-permanent use of AED training pads with embedded magnets in manikin (Offering magnetic stickers)




## SOFTWARE

- Download APP from Google Play
- Connect up to 6 manikins at the same time
- Editable guidelines
- Various Training mode : Game mode, Training mode, Evaluation mode
- Training for Chest Compression and Ventilation
  - Count Total Number and Correct Number of Compressions
  - Count Total volume and number of ventilation
  - Compression Depth and Speed
  - Incomplete recoil / Over-compression
- Save training result and score ; Export to Excel format
- Statistics & Comparison training data



\* Group Training



\* Result

## COMPONENTS

- Adult Torso : 1EA
- Tablet : 1EA
- Storage case (use for matress) : 1EA
- Power Adaptor : 1EA
- Magnetic stickers : 1SET (2EA)

## SPECIFICATION

- Dimensions : 65 x 33 x 20cm
- Weight : 4.2kg

