

THE PIONEER OF MEDICAL  
EDUCATIONAL SIMULATOR +

# Airway Management Training Model

BT-CSIE-T

The lifelike airway management training model

NURSING MEDICINE MILITARY EMS HOSPITAL



# Airway Management Training Model

## FEATURES

- Lifelike adult torso manikin
- Realistic structure and durable silicon skin
- Head tilt-Chin lift, Jaw thrust maneuver, Sniffing Position
- Insert NPA, OPA
- Insert ETT, LMA/i-gel, Combi-tube
- Intubation with a laryngoscope
- Ventilation with Bag valve mask
- Visible Chest rising / Gastric inflation during ventilation
- Check the wrong intubation with visible chest rising or gastric inflation
  - right mainstem intubation / esophageal intubation
- Breathing sound with stethoscope during ventilation



## COMPONENTS

- Adult Torso : 1EA
- Storage case : 1EA
- User manual : 1EA

## SPECIFICATION

- Dimensions : 70 x 30 x 24cm
- Weight : 8.5kg