

COMMODITY SPECIFICATION

Product	Baby First Aid Model	Model	BT-SEEM2-Baby
Manufacturer	BT Inc.	Qty	1 SET



A. Features

■ Software feedback via Bluetooth connection

- Supports Training mode / Evaluation mode
- Easily update guidelines and evaluation criteria
- Enables training for chest compressions, airway obstruction, and ventilation
- CPR skills can be trained and assessed using the software
- Displays Compression depth, rate, hand-position, ventilation volume and etc. on the S/W
- Simulated airway obstruction training supported via software
- Results can be reviewed, saved, and printed
- Generate statistics and comparative analysis based on assessment data
- Supports group evaluation by connecting up to 6 manikins to one tablet

COMMODITY SPECIFICATION

■ Hardware

- Visual chest rising while ventilation
- Embedded magnets in the manikin body, offering 2 magnetic stickers for AED pads
 - Use the AED trainer semi-permanently and without any damage on manikin's skin
- Replaceable face skin and lung bags

B. Components

- Baby full body manikin : 1 EA
- Software : 1 EA
- Lung bag, Airway tube : 1 EA/each
- Rubber steel sheets for AED pads : 2EA
- Power Adapter (12VDC 1.5A) : 1 EA
- Storage case : 1EA

-

C. Weight & Size

- Baby full body manikin : 3.4kg (350×670×160 mm)
- Storage case : 2.2kg (730 x 240 x 390 mm)
- Total package : Approximately 7 kg